Instructions for Use of Rotary Dental Instruments

Language	Document title	Page
English – EN	Instructions for use of rotary dental instruments	2
Deutsch – DE	Anwendungshinweise für rotierende Dentalinstrumente	6

BUSCH &CO. GmbH & CO. KG

Unterkaltenbach 17–27
51766 Engelskirchen
Germany

Phone: +49 2263 86-0
Email: mail@busch.eu
Website: www.busch.eu

CE 0124

Scope

These instructions for use apply to rotary dental instruments by BUSCH & CO. GmbH & Co. KG. They are intended exclusively for professional users in dental practices and dental clinics.

Intended Purpose, Indications, User, Patient Population

Intended Purpose

Rotary dental instruments

- Processing of dental hard and soft tissue, bone and dental restorative materials.

Polishers

- Smoothing and polishing dental hard tissue and dental restorative materials.

Indications

Steel and Carbide Instruments

- Processing (e.g. (pre-)grinding, finishing, removal, smoothing, correction and trepanation) of dental hard & soft tissue, bone and dental restorative materials.

Diamond Instruments

- Processing (e.g. (pre-)grinding, finishing, removal, smoothing, correction and trepanation) of dental hard & soft tissue, bone and dental restorative materials.

Depending on the specific instrument series, the indications also include:

- ZIRAMANT: Processing dental zirconia and ceramics
- Micro Prep: Minimally invasive preparation of hard-to-reach areas.
- Dia PERIO-PRO: Cleaning and smoothing of root surfaces.

Ceramic Abrasives

- Grinding of dental hard tissue and dental restoration materials.

Polishers and Prophylaxis Brushes

- Polishing for dentures of all restorative materials
- Polishing after placing a composite filling
- Polishing of precious and non-precious metal alloys, amalgam
- Polishing with paste for prophylaxis and shine-polishing of the enamel
- Brush for polishing of restorations made of ceramic, composite and metal
- Polishing of dental ceramics and silicate ceramics
- Polishing of dental restorations made of zirconia

Intended Patient Population

Dental burs are used for patients with damaged oral hard and soft tissue. Rotary dental instruments are used for prevention, prophylaxis, treatment or rehabilitation of the oral hard and soft tissue and dental restorations. Their main effect is mainly mechanical in nature. The intended target group consists of dental patients in dental offices and clinics.

Intended User

Dental rotary instruments are used in dental practices or clinics. The instruments are connected to a turbine or handpiece (e.g. straight or right-angle handpiece) and used by dentists or other professionals who are familiar with the safe handling of these instruments due to their training and experience.

Compatibility

Only transmission instruments suitable for rotary dental instruments of the respective shank type may be used. Determine the shank type from the product labelling. The shanks comply with EN ISO 1797. Observe the instructions for use of the turbine or handpiece used.

Contraindications

There are no known contraindications in connection with the devices.

Restrictions on Use and Pre-use Inspection

- Do not use for applications outside the stated dental intended purpose and indications.
- Observe product-specific material information. Where a relevant material intolerance or allergy is known or suspected, assess product suitability before use.
- Use only technically and hygienically sound, serviced and cleaned turbines, handpieces or contra-angle handpieces that are suitable for the intended purpose and compatible with the shank.
- Do not use an instrument if its shank, working part, coating, cutting edges or bond is damaged, corroded, deformed, worn or contaminated, or if the instrument cannot be clearly identified.
- Do not use the product if its packaging is damaged in a manner that may have affected product identity, integrity or cleanliness.
- Before use, perform a visual inspection for damage, wear and residual contamination. Dispose of heavily worn or damaged instruments. Re-clean contaminated reusable instruments according to the current reprocessing instructions available at www.busch.eu/eifu.

Warnings

Supplied Non-sterile and Single-use Products

WARNING: Reprocess reusable instruments before first use and immediately after each use in accordance with the current BUSCH reprocessing instructions for rotary instruments.

WARNING: Disposable products, which are identified by the symbol (⧖), may be processed only once before they are used for the first time and may not be reused. Reprocessing or reuse of single-use products may result in infection of the patient and/or user.

Maximum Permissible Speed

WARNING: Do not exceed the maximum permissible rotational speed indicated on the product label.

Turbines

WARNING: Use turbines only for the preparation of dental hard tissue. Do not use turbines for surgical procedures or crown cutting.

Application Force

WARNING: Bring the instrument to working speed before contact with the treatment site. Avoid tilting, levering, blocking and excessive application force. Apply only low pressure appropriate to the instrument and indication. Do not compensate for reduced performance by increasing the application force.

Cooling and Suction

WARNING: Ensure sufficient cooling appropriate to the instrument and indication. Stop use immediately if cooling is insufficient or excessive heat develops to avoid tissue damage.

Abnormalities During Use

WARNING: Stop use immediately in the event of unusual noises, vibration, run-out, unusual heating or an unexpected loss of performance. Inspect the instrument. Do not reuse a damaged or worn instrument; discard it.

Diamond Instruments

WARNING: The coating of diamond instruments contains nickel. Where nickel intolerance or allergy is known or suspected, assess product suitability before use.

Carbide Instruments

WARNING: Carbide instruments contain cobalt as a binder metal. Observe the applicable hazardous-substance information on the product labelling. Assess product suitability before use.

Side Effects and Complications

There are no known side effects or complications associated with the intended use of the devices.

Use

1. Select the instrument suitable for the intended application and verify its identity.
2. Inspect the instrument for contamination, wear, corrosion, deformation, damage, cracks, chipping, loosening, coating loss and other applicable discard criteria.

3. Securely clamp the instrument in the turbine or handpiece and verify secure seating and concentric running. The instrument must be clamped as deep as possible in the turbine or handpiece.
4. Select a suitable working speed below the maximum permissible speed indicated on the product label.
5. Bring the instrument to working speed before contact with material or tissue.
6. Guide the instrument in a controlled manner without tilting, levering or blocking. Apply only low application force (approx. 0.3 – 2N) appropriate to the instrument and indication.
7. Maintain cooling and suction appropriate to the indication. Where water cooling is required, provide at least 50 ml/min. Stop immediately if cooling is insufficient or excessive heat develops.
8. Do not compensate for loss of performance by increasing the application force. Inspect the instrument and replace it where necessary.
9. Allow the instrument to come to a complete stop before setting it down.
10. Transfer reusable instruments safely for reprocessing. Dispose of single-use products safely after the procedure.

Reprocessing

Reusable instruments must be reprocessed before first use and immediately after each use in accordance with the current BUSCH reprocessing instructions for rotary instruments available at www.busch.eu/eifu.

The specified procedures and parameters for cleaning, disinfection, drying, inspection, packaging and sterilisation must be observed. The detailed validated process parameters are provided in the current BUSCH reprocessing instructions.

Service Life and Reuse

Service life is determined by use, reprocessing, wear and the condition of the instrument.

For reusable instruments, reprocessing has been validated for at least 20 cycles. The validated number of reprocessing cycles does not constitute a general maximum number of permitted reuses. Reuse is permitted only after successful reprocessing and inspection and while the instrument continues to meet all specified acceptance criteria.

Reusable polishers are an exception: a maximum of 10 reprocessing cycles is permitted, provided that all inspection and discard criteria continue to be met.

Single-use products (⊗) may be processed only once before first use, as specified under “Supplied Non-sterile and Single-use Products”. After use, they must not be reprocessed or reused.

Discard an instrument no later than when any of the following applies:

- fracture, crack, chipping or loosening;
- bending, run-out, imbalance or abnormal vibration;
- corrosion or non-removable discolouration or contamination;
- blunt cutting edges, extensive grit loss, coating detachment or significant loss of performance;
- damaged shank or a shank that appears to be outside specification;
- brittle, cracked, deformed or loose polishing or abrasive bodies.

Storage and Transport

Store and transport instruments in a dry condition and protect them from dust, moisture, contamination and mechanical damage. Store processed instruments protected against recontamination.

Maintenance

No maintenance is required beyond the specified inspection and reprocessing measures.

Disposal

The instruments can be disposed of in normal household/practice waste. Dispose of packaging components in accordance with applicable national and local regulations.

Reporting Serious Incidents

Any serious incident that occurs in connection with the product must be reported to BUSCH & CO. GmbH & Co. KG and to the competent authority of the country in which the user and/or patient is established.

Please use the complaint form available in the download section at www.busch.eu/downloads and send it to mail@busch.eu.

Electronic Provision and Paper Copy

The instructions for use applicable to the respective product must be accessed at www.busch.eu/eifu. Before use, verify that the displayed version is assigned to the product and is identified as the currently valid version.

The electronic instructions for use are available for download and offline use. A free paper copy can be requested at mail@busch.eu and will be provided within the applicable regulatory period. Previous approved versions are archived in a controlled manner and clearly identified as not current.

Symbol Glossary

				
Manufacturer	Date of Manufacture	Importer	Distributor	Conformity marking
				
Batch Code	Catalogue Number	Unique Device Identifier	Medical Device	Consult electronic instructions for use
				
Do not re-use	Contains hazardous substances	Maximum permitted rotation speed per minute [min-1]	Authorized representative in the United Kingdom	Authorized representative in Switzerland

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