



Pusch®

There is no substitute for quality

# SIDE Grip



The SIDE GRIP a product out of our instrument range with passive safety.

The intelligent construction makes the new SIDE GRIP series to essential instruments for the selective abrasive finishing on infected thick nails or on critical parts of the skin.

SIDE GRIP instruments make treatments safer in edge areas preventing unintentional contacts.

New!

New!



805S



841S



840S





With our instrument range with passive safety, we offer pedicurists and podiatrists instruments with a special safety relevant construction for the practice.

A primary aim is to perform every treatment without any possible risk especially by treatments on diabetics. Our instrument range with passive safety supports the user by doing so.

Foremost by changes on infected nails and on certain skin parts the individual SIDE GRIP form can pinpoint individually the areas to be treated or rather to be abraded. Unintentional sliding contacts are avoided by the smooth frontal area.

We offer besides the approved cylindrical and conical basic form also the small dimensioned cylindrical form 841S with the option to do filigree treatments. The diamond coated conical SIDE-surface of the form 805S offers the advantage e.g. to abrade cornification safer in the front and on the sides of the nail edges.

**Recommended maximum attainable speed: 5.000 - 10.000 rpm**

Generally we recommend a decrease of rotation in the dry technique and also by treatments on high-risk patients.



**coarse grit**

**Safety bevel**

BUSCH	6840LS
Ø ISO	065

**medium grit**

BUSCH	840S	841S	854S	854S	840LS	805S
Ø ISO	055	035	033	040	065	065

**fine grit**

BUSCH	8840S	8841S	8805S
Ø ISO	055	033	065

= Package unit



**BUSCH & CO.** GmbH & Co. KG

Postfach 1152 • 51751 Engelskirchen • Unterkaltenbach 17-27 • 51766 Engelskirchen • GERMANY  
 Telefon +49 (0) 2263 860 • Telefax +49 (0) 2263 20741 • mail@busch.eu • www.busch.eu