





Hartmetall • Chirurgie-Instrumente, Knochenfräser
 Carbide • Surgical instruments, bone cutters
 Carbure • Instruments chirurgicaux, fraises à os
 Carburo • Fresas para la cirugía, fresas para huesos



161



FG	016
ISO 500 314 408 295...	
L mm	9,0



162



FG	016
ISO 500 314 408 297...	
H/HP/PM	016
ISO 500 104 408 297...	
W/RA/CA	016
ISO 500 204 408 297...	
W Ig/RA Ig/CA Ig	016
ISO 500 205 408 297...	
L mm	9,0



166



H/HP/PM	021
ISO 500 104 409 297...	
W/RA/CA	021
ISO 500 204 409 297...	
W Ig/RA Ig/CA Ig	021
ISO 500 205 409 297...	
L mm	10,0

Knochenfräser mit A-Verzahnung für schonendes Zerspanen und feinere Oberflächen

bone cutters with A-cut for a gentle chip removal and finer surface

fraises à os avec denture A pour un soigneux enlèvement des copeaux et une surface lisse


fresa para huesos con dentado A para arrancar cuidadosamente las virutas y finas superficies

Knochenfräser mit A-Verzahnung für minimalinvasives, schonendes Präparieren von Knochengewebe und Zahnhartsubstanz

bone cutters with A-cut for a minimal invasive and gentle preparation of bone tissue and hard tooth substance

fraises à os avec denture A pour un préparation minimalinvasif et soigneux du tissu osseux et de la substance dentaire dure

fresa para huesos con dentado A para preparación invasiva minimal y cuidadosamente del tejido de hueso y de la sustancia dura de dientes



162A




FG	016
ISO 500 314 408 298...	
H/HP/PM	016
ISO 500 104 408 298...	
W/RA/CA	016
ISO 500 204 408 298...	
W Ig/RA Ig/CA Ig	016
ISO 500 205 408 298...	
L mm	9,0



166A



H/HP/PM	021
ISO 500 104 409 298...	
W/RA/CA	021
ISO 500 204 409 298...	
W Ig/RA Ig/CA Ig	021
ISO 500 205 409 298...	
L mm	10,0



164LA



FG	012
ISO 500 314 415 298...	
L mm	6,5